

## **ENUMERATION OF BACTERIA, YEASTS AND MOLDS: STANDARD PLATE COUNT**

The standard plate count method consists of diluting a sample with sterile saline or phosphate buffer diluent until the bacteria are dilute enough to count accurately. That is, the final plates in the series should have between 30 and 300 colonies. Fewer than 30 colonies are not acceptable for statistical reasons (too few may not be representative of the sample), and more than 300 colonies on a plate are likely to produce colonies too close to each other to be distinguished as distinct colony-forming units (CFUs).

### **MATERIALS NEEDED**

- 1) 24-hour 10 ml nutrient broth culture of *Escherichia coli*
- 2) 4 sterile 99-ml saline bottles
- 3) 1ml and 5ml pipettes
- 4) 6 Petri plates
- 5) 6 agar pour tubes of nutrient agar (plate count agar)

### **PROCEDURE**

1. Using aseptic technique, the initial dilution is made by transferring 1 ml of *E. coli* sample to a 99ml sterile saline blank. This is a  $1/100$  or  $10^{-2}$  dilution.
2. Immediately after the  $10^{-2}$  dilution has been shaken, uncap it and aseptically transfer 1ml to a second 99 ml saline blank. This second blank represents a  $10^{-4}$  dilution of the original sample.
3. Shake the  $10^{-4}$  dilution vigorously and transfer 1ml to the third 99 ml blank. This third dilution represents a  $10^{-6}$  dilution of the original sample. Repeat the process once more to produce a  $10^{-8}$  dilution.
4. Shake the  $10^{-4}$  dilution again and aseptically transfer 1.0 ml to one Petri plate and 0.1 ml to another Petri plate. Do the same for the  $10^{-6}$  and the  $10^{-8}$  dilutions.
5. Remove one agar pour tube from the 48 to 50<sup>0</sup> C water bath. Carefully remove the cover from the  $10^{-4}$  Petri plate and aseptically pour the agar into it. The agar and sample are immediately mixed gently moving the plate in a figure-eight motion or a circular motion while it rests on the tabletop. Repeat this process for the remaining five plates.
6. After the pour plates have cooled and the agar has hardened, they are inverted and incubated at 25<sup>0</sup> C for 48 hours or 37<sup>0</sup> C for 24 hours.
7. At the end of the incubation period, select all of the Petri plates containing between 30 and 300 colonies. Plates with more than 300 colonies cannot be counted and are designated too many to count (TMTC). Plates with fewer than 30 colonies are designated too few to count (TFTC). Count the colonies on each plate. A Quebec colony counter is used to count the colonies in the culture plate.
8. Calculate the number of bacteria (CFU) per millilitre or gram of sample by dividing the number of colonies by the dilution factor multiplied by the amount of sample added to liquefied agar.

$$\text{Number of No. of bacteria} = \text{No of colonies (CFUs)} / \text{dilution} \times \text{amount plated}$$

## Direct microscopic count

This method requires a haemocytometer (Petroff-Hausser chamber)

### Batch Growth

- Quantity of cell concentration — determination of cell number density


The diagram illustrates the direct microscopic count method. It shows two main components: the Petroff-Hausser slide and the hemocytometer. The Petroff-Hausser slide consists of a cover glass and a chamber holding bacteria. The hemocytometer consists of a cover glass support and a counting chamber with a depth of 0.1 mm. The diagram shows the process of counting cells in a grid, with a red grid representing the Petroff-Hausser slide and a black grid representing the hemocytometer. The counting chamber is shown with a depth of 0.1 mm.

Petroff-Hausser slide


hemocytometer

Drug Delivery & Tissue Engineering lab.  
Sogang University

서강대학교  
SOGANG UNIVERSITY



Depth between counting chamber and cover slip is mostly 0.1 mm



**0.2 mm**

### **Large square dimensions:**

Area:  $1 \times 1 \text{ mm}^2$  Volume:  $1 \times 1 \times 0.1 \text{ mm}^3$  or  $(1 \times 1 \times 0.1)/1000 \text{ cm}^3$

Calculation:  $n \times 10^4$

### **Smaller square dimensions:**

Area:  $(1 \times 1)/(5 \times 5) \text{ mm}^2$  Volume:  $[(1 \times 1)/(5 \times 5)] \times 0.1 \text{ mm}^3$  or  $[(1 \times 1)/(5 \times 5)] \times 0.1/1000 \text{ cm}^3$

Calculation:  $n \times (25 \times 10^4)$

### **Smallest square dimensions:**

Area:  $(1 \times 1)/(20 \times 20) \text{ mm}^2$  Volume:  $[(1 \times 1)/(20 \times 20)] \times 0.1 \text{ mm}^3$  or  $[(1 \times 1)/(20 \times 20)] \times 0.1/1000 \text{ cm}^3$

Calculation:  $n \times (400 \times 10^4)$  or  $n \times (4 \times 10^6)$

### **Procedure**

Methylene blue is added to the sample of bacterial suspension.

Petroff Hausser chamber's grid is covered with a cover slip and a drop of sample is placed in the chamber with a Pasteur pipette.

Sample observed under microscope.


Count the number of cells in five smaller squares.

Calculate the average number of cells for five squares.

Estimate the number of cells (cfu) for one ml by the formula

**Cfu/ml = average number of cells in smaller square X (25 X 10<sup>4</sup>) X dilution factor**

## Direct microscopic count



In direct microscopic counting: 1) dead cells are not distinguished from living cells; 2) small cells are difficult to see under the microscope; 3) precision is difficult to achieve; 4) we need a phase contrast microscope; 5) not a good method for cell suspensions of low density.